



PetGuard Pro

Hardware User Manual

Smart Sensor Pod & Harness — A Guide for Pet Owners

Prepared by the PetGuard Pro Team
(e21-3yp-PetGuard-Pro / Hardware)

Contents

1. Introduction
2. What's in the Box
3. Getting to Know Your Sensor Pod
4. Attaching the Sensor Pod to the Harness
5. Fitting the Harness on Your Pet
6. Powering the Device On and Off
7. Charging the Device
8. The Temperature Sensor — Important Care Notes
9. Safety Precautions & Care
10. Connecting to the PetGuard Pro App
11. Cleaning & Maintenance
12. Troubleshooting
13. Technical Specifications
14. Support & Contact

1. Introduction

Thank you for choosing PetGuard Pro. This manual covers the physical hardware that ships with your PetGuard Pro system — the Sensor Pod, the harness, and the charging accessories — and shows you how to set them up, wear them safely, and keep them working well.

For instructions on using the PetGuard Pro mobile app (creating an account, viewing live vitals, setting safe zones, and more), please refer to the separate PetGuard Pro Mobile Application User Manual.

Your PetGuard Pro hardware lets you:

- Track your pet's live GPS location and receive safe-zone alerts
- Monitor body temperature in real time
- Track daily activity — steps, rest, and sudden impacts such as falls
- Send all of this data automatically to the PetGuard Pro app over the mobile network

2. What's in the Box

Your PetGuard Pro package includes:

- 1 × PetGuard Pro Sensor Pod (the white enclosure with the on/off switch and temperature sensor)
- 1 × Adjustable pet harness with a built-in Velcro mounting panel
- 1 × USB charging cable
- 1 × Quick start card
- This Hardware User Manual

If any of these items are missing or appear damaged, please contact PetGuard Pro Support before using the device — see Section 14.

3. Getting to Know Your Sensor Pod

The Sensor Pod is the small white unit that houses all of PetGuard Pro's electronics — GPS, motion sensor, temperature sensor, cellular module, and battery. It clips onto the harness using a Velcro strap so it can be removed for charging or cleaning.

3.1 Front Panel — On/Off Switch



The red rocker switch on the side of the Sensor Pod — press “I” for on, “O” for off.

3.2 Underside — Temperature Sensor



The small circular sensor in the centre of the underside is the body-temperature sensor.

3.3 End Panel — Charging Port



The USB charging port and power switch are located on the end of the Sensor Pod.

3.4 Charging Status LED

Next to the switch is a small LED indicator that shows the battery's charging status whenever the Sensor Pod is connected to a charger:



Red LED — the battery is currently charging.



Blue LED — the battery is fully charged.

- Red light: charging is in progress. Leave the device connected to the charger.
- Blue light: charging is complete. It's safe to unplug the Sensor Pod.

4. Attaching the Sensor Pod to the Harness

The harness has a dedicated Velcro mounting panel on the back pad. Follow these steps to attach the Sensor Pod:

1. Lay the harness flat with the back pad (the wide padded section) facing up.
2. Locate the two Velcro strips on the back pad.
3. Place the Sensor Pod on the Velcro panel with the temperature sensor (the small round window) facing down, so it will sit against your pet's skin/fur when worn.
4. Press down firmly along both Velcro strips until the pod is fully seated and does not wobble.
5. Give the pod a gentle tug to confirm it is securely fastened before fitting the harness on your pet.



Sensor Pod and Velcro panel before attachment.



Sensor Pod correctly attached to the harness.

Tip

Always attach the Sensor Pod before putting the harness on your pet — it's much easier to seat the Velcro properly on a flat surface.

5. Fitting the Harness on Your Pet

1. Loosen all adjustable straps before putting the harness on.
2. Slip the harness over your pet's head so the padded chest section sits centrally on the chest.
3. Bring the belly strap underneath and clip it into the side buckle.
4. Tighten the shoulder and belly straps evenly until the harness is snug but not tight — you should be able to fit two fingers flat between the strap and your pet's body.
5. Check that the Sensor Pod sits flat against your pet's side or back, with the temperature sensor window making gentle contact with fur/skin and no strap covering it.
6. Let your pet walk around for a minute and re-check the fit — straps can shift slightly on first wear.



Harness top view, showing the adjustable straps and D-ring.

Fit Check

A properly fitted harness should not slide sideways or ride up toward the neck when your pet walks. Recheck the fit weekly, and always after a growth spurt for puppies.

6. Powering the Device On and Off

The Sensor Pod has a single rocker switch on its side panel.

- To turn ON: press the side marked “I”.
- To turn OFF: press the side marked “O”.

Turn the device on a minute or two before putting the harness on your pet, so the GPS and cellular module have time to connect. There is no separate power-on chime; check the PetGuard Pro app's Collar Status card to confirm the device is Online.

Battery Saving

Always switch the Sensor Pod off when it is not being worn, such as overnight in storage or while charging is not required, to preserve battery life for your next walk.

7. Charging the Device

1. Switch the Sensor Pod off before charging.
2. Locate the USB charging port on the end panel of the Sensor Pod, next to the power switch.
3. Connect the supplied USB cable to the port and to a standard USB power source (5V).
4. While charging, the LED next to the switch glows red. Charging typically takes 1.5–2 hours from empty.
5. When the LED turns blue, the battery is fully charged. Unplug the cable and store the Sensor Pod somewhere dry.

Warning

Only use the supplied cable or a standard 5V USB charger. Do not charge the device while it is being worn by your pet. Do not leave the device charging unattended for long periods.

8. The Temperature Sensor — Important Care Notes

The small circular window on the underside of the Sensor Pod is the body-temperature sensor. It must be able to sense your pet's body temperature accurately, so please follow these guidelines carefully:

Do Not Touch or Cover

Never press, poke, or place tape or stickers over the temperature sensor window. Doing so can damage the sensor or produce inaccurate readings.

- The sensor must stay in gentle contact with your pet's fur/skin — this is why the Sensor Pod is mounted with the sensor window facing inward, toward the body.
- Do not fit the harness with the Sensor Pod upside down or sideways, as this will lift the sensor away from your pet and cause inaccurate or missing temperature readings.
- Do not place any padding, clothing, or accessories between the Sensor Pod and your pet's body.
- If the app shows an unusually high or low temperature reading, first check that the sensor is making proper contact before assuming a health issue — then consult the Troubleshooting section or a veterinarian if the reading persists.

9. Safety Precautions & Care

9.1 General Safety

- PetGuard Pro is designed to complement, not replace, direct supervision of your pet.
- Do not leave the harness on an unsupervised pet for extended periods; remove it during rest, crate time, and sleep.
- Check the harness fit regularly, especially for growing puppies, to avoid chafing.
- Keep the charging cable and small accessories away from pets and children when not in use.

9.2 Sensor Pod Care

- The Sensor Pod is splash-resistant but is not designed for submersion — avoid swimming, bathing, or heavy rain exposure while it is attached.
- Avoid dropping the Sensor Pod or exposing it to extreme heat (such as a parked car) or direct sunlight for long periods.
- Remove the Sensor Pod before washing the harness.

9.3 Battery Safety

- Use only the supplied cable and a standard 5V charger.
- Stop using the device and contact support if you notice swelling, unusual heat, or a damaged casing.
- Do not attempt to open the Sensor Pod casing — there are no user-serviceable parts inside.

10. Connecting to the PetGuard Pro App

The Sensor Pod does not need to be manually paired — once your PetGuard Pro account has been approved and a Sensor Pod has been assigned to your pet by an administrator, the device will connect automatically over the cellular network whenever it is switched on and within signal range.

1. Switch the Sensor Pod on.
2. Open the PetGuard Pro app and go to the Home tab.
3. Wait up to 1–2 minutes for the Collar Status card to change from Offline to Online.
4. Once online, live vitals and GPS coordinates will begin updating automatically.

For account setup and app navigation, refer to the PetGuard Pro Mobile Application User Manual, Sections 2–4.

11. Cleaning & Maintenance

- Remove the Sensor Pod before cleaning the harness.
- Wipe the Sensor Pod casing with a soft, slightly damp cloth. Do not submerge it or use solvents.
- Hand-wash the harness fabric in cold water with mild soap; air-dry fully before reattaching the Sensor Pod.
- Gently clean the temperature sensor window with a dry, soft cloth if it appears dirty — avoid scratching it.
- Inspect straps, buckles, and Velcro periodically for wear and replace the harness if any part is frayed or the buckle no longer clicks securely.

12. Troubleshooting

Issue	What to try
Device won't turn on	Charge the Sensor Pod for at least 30 minutes, then try the switch again. Confirm the switch is fully pressed to "I".
harness shows Offline in the app	Confirm the Sensor Pod is switched on and charged. Check that you are in an area with mobile network coverage. Restart the Sensor Pod (switch off, wait 5 seconds, switch on).
No GPS location updates	Make sure the Sensor Pod has a reasonably unobstructed view of the sky. Allow up to 2 minutes for GPS to acquire a signal outdoors from a cold start.
Temperature reading looks wrong or shows "--"	Check the harness fit — the sensor window must sit flat against your pet's body with no strap or padding in between. Refit the harness and wait a few minutes for a fresh reading.
Heart Rate always shows "--"	This is expected — the current Sensor Pod does not include a heart-rate sensor (see Section 1).
Sensor Pod feels warm while charging	This is normal during charging. If it becomes hot to the touch or swells, disconnect immediately and contact support.

LED stays red for a very long time	Confirm the cable and power source are working (try a different USB adapter/cable). A full charge from empty normally takes 1.5–2 hours.
No LED light at all when plugged in	Check the cable is fully inserted at both ends and the power source is switched on. Try a different USB adapter or cable.
Sensor Pod won't stay attached to the harness	Check that both Velcro surfaces are clean and free of pet hair or debris. Press down firmly along the full length of the strap when reattaching.

For anything not covered here, please contact PetGuard Pro Support — see Section 14.

13. Technical Specifications

Feature	Details
Sensors	GPS module, 3-axis motion sensor (activity & impact detection), body-temperature sensor
Connectivity	Cellular (2G/GPRS) data connection to the PetGuard Pro cloud
Battery	Rechargeable Li-Po battery, USB charging, approx. 1.5–2 hour full charge
Typical battery life	Approx. 1 full day of normal, active use per charge
Enclosure	Splash-resistant 3D-printed casing; not designed for submersion
Mounting	Velcro-panel attachment to the included adjustable harness
Operating conditions	Avoid prolonged exposure to direct sunlight, extreme heat, or heavy rain
Heart rate sensing	Not included in the current hardware revision
Charging status LED	Red = charging in progress; Blue = fully charged

14. Support & Contact

If you experience a hardware issue that isn't resolved by the Troubleshooting section, please reach out through the Help & Support screen in the PetGuard Pro app, or contact the PetGuard Pro team directly with your Pet ID and a description of the issue.

Please keep this manual for future reference and share it with anyone else who cares for your pet.

 *Thank you for trusting PetGuard Pro with your pet's safety.* 